

# FlowLab

Software for CITREX and FlowAnalyser

analyser  
the art of measuring

## Reporting and documentation software

View your flow, pressure and volume measurements comfortably on your computer monitor while simultaneously saving respiratory parameters. FlowLab displays your measurement results graphically in real-time curves or numerically. Trending reports of up to 100 hours can be generated to verify long-term ventilator functionality. To facilitate your device management, test reports can also be created and saved electronically or printed.

« Software package for real-time curves, 100h trending and customized reporting. »»



## Test Report

- Test report printouts with one click
- User defined configuration
- Logo insertion available
- FlowAnalyser data automatically retrieved
- Various input options for each tested ventilator
- Unique control number for each report

## Panels

- Displayed in time-relation or as a loop
- Various cursors measure the curves
- Unique trigger used to display real-time curves in Single Shot, Norm or Auto mode
- User defined layout and colours
- Option of setting a title, printing and saving
- Simultaneous display of up to 6 different curves

## Trending

- Up to 100 hours of data logging
- User defined trending interval
- Up to 10 values simultaneously
- Data export to Excel, etc.
- User defined layout and colours
- Option of setting a title, printing and saving
- Automatic file size calculation

## Numerics

- All measurements displayed on one page or combined with panels
- User defined layout and colours
- Statistical data including mean, max and min for each value
- Input of target value with tolerances
- Automatic verification of Pass/Fail
- Up to 20 values can be displayed simultaneously

## Languages

- Available in 15 languages
- All languages included – free of charge
- Languages from every continent

